


QTinno™ Overview



 NewCardio 3-Dimensional software for measuring and adjudicating all timing intervals and voltage/amplitude parameters

 Key features include:

- All measurements made on the heart **Vector Magnitude lead**
- **Novel, proprietary algorithms** for accurate and precise fiducial point placement
- **Full automation** with intelligent self-validation algorithms to ensure reliability
- **Confidence Factor™** identifies and flags problematic ECGs for possible human overread
- **Analysis of all readable P-QRS-T complexes** in each ECG or Holter recording
- **Continuous quantitative assessment** of ventricular repolarization
- **No change in ECG practice**; accepts all files from standard ECG Equipment

Key Features of QTinno Engine


- Use of Vector Magnitude Virtual Lead
 - Produces single “virtual” ECG lead with full electrical information
 - from *all* 12 leads and *all complexes in an ECG*
 - improves signal-to-noise ratio over standard 12-lead ECG
- Proprietary, iterative curve fitting algorithms optimize fiducial point accuracy and precision
 - Least-square fit of 3rd-order polynomial to VM data
 - Multiple iterations done to optimize fitting
 - Minimum of fitted function represents T_{end}
 - Approach is noise-resistant and baseline-independent



Role of Confidence Factor™

- Internal assessment of quality and reliability is a mandatory feature of automated ECG analysis - otherwise, user must choose between 2 undesirable options: (1) accept automated results on a “black box” basis, or (2) overread all results
- Confidence Factor assigns each ECG a CF score of 0 to 100 based on:
 - ECG signal quality, including excessive noise, missing leads, etc
 - Measurement quality, including goodness of fit, heart rate variability, etc
 - Currently, $CF < 70$ flags an ECG for possible overread (about 1-5% of early-phase studies)

*Mean Std Dev of Difference
Between QTinno and Core Lab,
by CF Category*





QTinno™ Performance in Clinical Validation Studies

Comparison Between Methods (QTinno – Core Lab)



Study	CPU time	Diff in Mean Δ QTcF	STD of Diff, Raw QTcF	Categorical Analysis	ECGs flagged for overread by QTinno
NCE1 (n=1963)	<10 min	<0.1 ms	<7 ms	Highly concordant	<2 %
NCE2 (n=7039)	<30 min	<1 ms	<8 ms	Highly concordant	<2 %
NCE3 (n=2523)	<15 min	<1 ms	<10 ms	Highly concordant	<10% (high-noise dataset extracted from Holters)
TQT1 (n=10,880)	< 30 min	<0.1 ms	<8 ms	Highly concordant	<5%
TQT2 (n=17,506)	< 60 min	<1 ms	<8 ms	More outliers reported by Core Lab than by QTinno	<1%
TQT3 (n=29,547)	< 120 min	Analysis in Progress			
TQT4 (24 Hr Holters , 72 pts)	< 60 min per Holter	Analysis in Progress			

Comparison of 2 TQT Studies - Methods

TQT1 Study


-  Pharma-sponsored dataset from parallel-group TQTS; 10,880 ECGs from moxi and placebo arms
-  QTinno fully automated results generated in blind fashion, then analyzed and compared by Pharma physician-scientists to manual Core Lab measurements

TQT2 Study

-  Unblinded CSRC dataset from crossover TQTS; 17,506 ECGs from moxi and placebo arms
-  QTinno fully automated results internally compared to Core Lab results


Comparison of 2 TQT Studies – ddQTcF by Hour

Mean ddQTcF



TQT1 Study

- Both Core Lab and QTinno show typical moxi ddQTcF curve
- Less than 1 ms difference in point estimates at 7 of 11 time points
- 90% CIs for QTinno and Core Lab virtually identical in width




TQT2 Study

- Both Core Lab and QTinno show typical moxi ddQTcF curve
- Less than 1 ms difference in point estimates at 7 of 12 time points
- 90% CIs consistently narrower for QTinno relative to Core Lab

Comparison of 2 TQT Studies – Method Differences in dQTcF


Bland-Altman Analysis
Core Lab vs. QTinno Automated

TQT1



Mean difference = **-0.1 ms**
Std Dev = 8.1 ms

TQT2



Mean difference = **0.0 ms**
Std Dev = 9.8 ms

Comparison of 2 TQT Studies: Categorical (Outlier) Analysis for dQT/dQTcF

TQT1 Study

Number of Subjects by Category

	<30 ms	30-60 ms	60-90 ms	>90 ms	Total
--	--------	----------	----------	--------	-------

dQT		Number of Subjects by Category				Total
		<30 ms	30-60 ms	60-90 ms	>90 ms	
Plac	Core Lab	437	27	5	0	469
	QTinno	434	28	6	1	469
Moxi	Core Lab	458	45	7	0	510
	QTinno	452	48	9	1	510

dQTcF		Number of Subjects by Category				Total
		<30 ms	30-60 ms	60-90 ms	>90 ms	
Plac	Core Lab	466	3	0	0	469
	QTinno	466	3	0	0	469
Moxi	Core Lab	495	15	0	0	510
	QTinno	489	21	0	0	510

TQT2 Study

Number of Subjects by Category

	<30 ms	30-60 ms	>60 ms	Total
--	--------	----------	--------	-------

	Number of Subjects by Category			Total
	<30 ms	30-60 ms	>60 ms	
Core Lab	1396	40	1	1437
QTinno	1403	24	0	1437
Core Lab	1385	76	0	1461
QTinno	1425	36	0	1461

	Number of Subjects by Category			Total
	<30 ms	30-60 ms	>60 ms	
Core Lab	1434	3	0	1437
QTinno	1437	0	0	1437
Core Lab	1434	27	0	1461
QTinno	1459	2	0	1461


📍 Shows substantial agreement between Core Lab and QTinno in both categorical analysis and distribution of dQT/dQTcF, for both Placebo and Moxifloxacin groups

File Filters View

12 Leads	QT View	Global	10 sec	QT (3C)	QT (All)	PR	QRS	Unlock	< Prev.	Next >	Reject Tracing	Accept Read...			
I	II	III	aVR	aVL	aVF	V1	V2	V3	V4	V5	V6	X	Y	Z	V

QT Client 2 - Java - version 1.075i - NewCardio Proprietary - use only under signed CDA!

QT1a: 371 QT2a: 364 QT3a: 362 CT: 3.064 RR: 0.872 sec CFT: 100 CFR: 11 CFF: 1 CF: 92 RRVar=N QT: 365 ms QTcf: 362 ms ORR: no ADJ: acc by: n/a GS: 321 (ms) PR: 162 PRCF: 60 QRS: 70 QRSCF: 93



Residual Variability Obtained from Fitting a Mixed Model to dQTcF by Method


Study	Method	Intrinsic Variability (ms)	
		Between-Subject	Within-Subject
TQT1	Core Lab	7.4	8.8
	QTinno	7.9	8.1
TQT2	Core Lab	3.8	10.7
	QTinno	3.2	6.7

 Mixed Statistical Model with factors for day, time, age, gender, treatment and all interactions

Holter-Based TQT Study Design

- 📍 Moxi and Placebo arms from a recent Pharma-sponsored TQTS
- 📍 Continuous 12-Lead Holter recordings on 71 healthy volunteer subjects
- 📍 Crossover design; all subjects have 1 day on placebo and 1 day on moxi
- 📍 Baseline for both placebo and moxi days – average of hrs -1, -0.5, and 0 for all subjects
- 📍 Ten predetermined time points for 12-lead ECG extraction – hrs 1, 2, 3, 4, 5, 6, 8, 12, 23, 24
- 📍 Extraction performed within 5 min of predetermined time point, after 10 min of supine rest


Hourly ddQTcF by Measurement Method



📍 All 3 QTinno measurement methods produced acceptable and similar ddQTcF moxi curves, with All Good Beats (shown) having the narrowest CIs

📍 Results closely comparable to Core Lab results

“All Good Beats” vs QTc-Based Selection: Effect on Number of Beats Analyzed



- Using QTc-based 10 min Screening Window, 30 sec Measurement Window, CF>70 cutoff results in an average of 31 beats analyzed per subject per datapoint
- Using “All Good Beats”, 10 min window, CF>70 cutoff results in 466 beats analyzed per subject per datapoint
- This may be an advantage for AGB that contributes to substantially lower 90% CIs than QTc-based selection

Residual Variability Obtained from Fitting a Mixed Model to dQTcF by Method

Method	Intrinsic Variability (ms)	
	Between-Subject	Within-Subject
Core Lab	TBD	6.9
QTinno (RR-based selection)	3.9	6.2
QTinno (QTc-based selection)	4.3	6.0
QTinno (All Good Beats)	3.5	5.4

 Mixed Statistical Model with factors for day, time, age, gender, treatment and all interactions

Summary and Conclusions

- 🌀 NewCardio automated analysis yields acceptable ddQTcF moxi curves for all measurement approaches
 - 🌀 “All Good Beats” appears to yield better data than either QTc-based or RR-based screening with 30 or 60 sec measurement windows
 - 🌀 Cutoff of $CF > 70$ is superior to $CF > 0$ cutoff, for QTc-based or RR-based screening, but CF cutoffs were approximately equal for All Good Beats
 - 🌀 For All Good Beats, larger measurement windows appear superior to smaller, and leftward-shifted measurement windows appear superior to rightward-shifted windows
- 🌀 NewCardio automated analysis performs well on measurement of non-QT IDMs, such as RR, PR, and QRS intervals