



In Vitro Diagnostic Testing for Direct Oral Anticoagulants- Premarket Review

Cardiac Safety Research Consortium

December 3, 2015

Lea Carrington

Division of Immunology and Hematology Devices
Office of In Vitro Diagnostics and Radiological Health
Center for Devices and Radiological Health

Objectives

- Overview of the regulatory requirements for clearance of IVDs intended for coagulation testing in patients treated with DOACs
- Identify clinical circumstances/conditions where testing of DOACs anticoagulant activity or concentration would be relevant
- Understand limitations of assays used to monitor novel anticoagulants

FDA Regulated Uses of IVDs

- **Diagnosis** – Diagnose disease, identify, confirm or rule out a condition in symptomatic patients
- **Screening** – Intended use population includes individuals **without** signs or symptoms of disease
- **Monitoring, prognosis, prediction**

What kind of submission should be prepared?

Classification of IVD Devices

- Regulatory path determined using a risk-based approach
 - Class I: low likelihood of harm, most are 510(k) exempt
 - Class II: Premarket notification, 510(k)
 - Class III: High or unknown likelihood of harm or how to prevent harm is unknown – Premarket Approval (PMA)
- Classification (I, II, or III) depends on risk

Risk Depends on Intended Use

- Risk (and subsequently device classification and submission type) is inherently tied to the claimed Intended Use of a device
- Level of FDA review and type of studies requested generally depend on the IU claims



Terminology

Class	Premarket Submission	Success Metric	Action/Decision
III	PMA	Safety and Effectiveness	Approval
II	510(k)	Substantial Equivalence	Clearance
I	None (if exempt)		
II (De novo)	510(k)	Safety and Effectiveness	Clearance

FDA Regulatory Considerations

- Intended Use/Indications for Use - who should be tested?
- What clinical evidence is needed to support claims?
- Can effectiveness of test system be reliably determined from data/information provided?
- What are the analytical performance requirements?
- Is it necessary to restrict use of test system – settings/conditions?
- How is the device output to be clinically interpreted?

Review Elements of IVD Submissions

- **Intended Use: how and by whom the device is to be used**
 - Measurand (analyte, biological activity)
 - Testing population
 - Quantitative, Semi-quantitative or Qualitative
 - Specimen type(s)
 - Conditions for use/clinical setting/intended user
- **Indications for Use: for what and for whom the device is to be used**
 - Conditions or disease
 - Target population (e.g. age)
 - Time and frequency of use

Performance Characteristics

- **Analytical performance** – does the test device measure the analyte correctly? How reliably?
- **Clinical performance** – does the test result correlate with the expected clinical presentation? How reliably?
 - Appropriate study design, description?
 - Samples tested near cut-off or relevant medical decision points?
 - Informative in the context of intended use?
 - All sample matrices evaluated?
 - Pre-analytical steps included in evaluation?

Review Elements of IVD Submissions

- Performance characteristics: **Analytical**
 - Establish basic performance parameters
 - Precision/reproducibility, reference range/cut-off, analytical sensitivity & specificity
 - Acceptance criteria (established *a priori*)
 - Are the established limits reasonable to support validation of the proposed intended use?
 - What are the appropriate studies to evaluate performance as defined in the intended use?
 - Use of traceable reference materials and methods, if available; calibrators and controls

Review Elements of IVD Submissions

- Performance characteristics: **Clinical**
 - Study design should include the target population
 - Prospective vs. Retrospective samples
 - Matrix considerations, claimed sample type(s)
 - Number of sites
- Comparator
 - Predicate or Reference method
 - Clinical performance

Key Issues

- Manufacturers are developing devices to assess effect or concentration of Direct Oral Anticoagulants (DOACs)
- Acceptable result – numerical value/qualitative detection
- Monitoring is not routinely required
- Currently no cleared or approved devices

FDA's Goal: To ensure that we enable timely access to safe, effective, and high-quality medical devices

Regulatory Challenges

- Lack of good correlation with drug concentration; relation to clinical outcomes
- Adequacy of expected values, reference ranges to support regulatory assessment for clinical use
- Appropriate time interval required to assess drug effect in coagulation assays post dosing
- Availability of published data on therapeutic levels

Clinical Laboratory Testing

- Potential situations where laboratory measurement of DOAC drug effect may be useful:
 - Assess compliance
 - Hemorrhagic or thrombotic complications
 - Determination of overdose
 - Pre-operative, urgent surgery

Early Interactions with CDRH

- Utilize the Pre-submission process
 - Highly recommended, offers free FDA consultation
- Provide a clear, carefully crafted intended use
- Provide protocols with as much detail as possible, include specific questions
- Make sure proposed studies support IU claims
- Open communication with the Agency to discuss regulatory strategy and study design



Questions?

Lea.Carrington@fda.hhs.gov